

Osteomyelitis Toe

(revised 2/16)

Contrast dosage 0.1mmol/kg

Mark area(s) of interest(s). Do not depress capsule into skin, place superficially

If the area of concern is greater than the toe then please scan thru that area as well

If area of interest ambiguous call/page the resident or attending

If poor T2 fat saturation use STIR

For Toe narrow FOV to area

1. SCOUT

2. SHORT AXIS T1 TSE **2mm only thru toe/bone of interest**

3. SAG T1 TSE “

4. SAG STIR “

5. SHORT AXIS T2 FS/STIR “

6. LONG AXIS T2 FS FOV **max 14**

7. SHORT AXIS T1 FAT SAT **2mm only thru**

Inject Contrast

8. SAG T1 TSE FS **2mm only thru toe/bone of interest**

9. SHORT AXIS T1 TSE FS “

10 LONG AXIS T1 TSE FS **larger FOV max 14 and thicker slices ok**