

# **Athletic Pubalgia**

## **Center on symphysis pubis**

Use body array coil

Have patient empty bladder prior to scan

\*pulsation artifacts NOT thru pubic symph – swap phase and freq

Axial plane small fov: scan from above the acetabulum to the lesser trochanter

Sagittal plane small fov: scan perpendicular to coronal plane through both ischial tuberosities

Coronal plane small fov: prescribe plane with slices parallel to the symphysis pubis, scan through the ischial tuberosities and symphysis pubis

## **1. SCOUT**

## **LARGE FOV COVER ENTIRE PELVIS**

**2. COR PD HYBRID fs** - call to check

**3. COR SE T1**

## **SMALL FOV on pubic symphysis**

**4. SAG PD HYBRID fs**

**5. AX PD HYBRID fs**

**6. COR PD HYBRID fs**

**7. AX PD (high resolution) non fs**