

Forefoot (Non osteomyelitis)

(revised 12/ 2009)

For sagittal imaging use a long and short axis scout to localize. Orientate slices perpendicular to metatarsals and parallel to the long axis of the foot. Cover to the air medial to lateral

For long axis imaging cover from the medial to the lateral aspect of the foot

For short axis imaging cover from the toes to the tarsals

1. SCOUT

2. SAG TSE T1

3. SAG STIR

4. Short axis TSE T1 (perpendicular to metatarsals)

5. Short axis PD HYBRID FS

6. Long axis PD HYBRID

7. Long axis PD HYBRID FS