

KNEE 3T

1. For axial imaging cover from the insertion of the patellar tendon to 2cm above the patella or through any fluid collection.
 2. For sagittal imaging, angle perpendicular to femoral condyles
 3. For coronal imaging, cover from the neurovascular bundle posteriorly through the patella anteriorly with the slices parallel to the femoral condyles.
 4. Use the extremity coil
- Revised 06/2014

1.Scout

2.Ax PD Hybrid fs

3.Sag T2 fs

4.Sag PD Hybrid non fs

5.Cor PD Hybrid fs

6.Cor PD Hybrid non fs

7.Cor SE T1

OPTIONAL

Sag WATSf only if cartilage sequences requested by ordering physician, or requested by radiologist