

Liver

(revised 12/2009)

Body phase array coil.

Use power injector

Cover entire liver and pancreas

Contrast dosage 0.1mmol/kg

Max slice thickness 6mm

- 1. SCOUT BH**
- 2. AXIAL FRTSE T2 respiratory compensation**
- 3. AXIAL FRTSE T2 FS resp comp**
- 4. AXIAL T1 GRE OUT/IN PHASE**
- 5. AXIAL DIFFUSION (B500/B0)**
- 6. AXIAL LAVA pre**

Inject – Fluoro-prep

- 7. AXIAL POST GAD LAVA BH dynamic**
- 8. CORONAL POST GAD LAVA BH (2 minute)**
- 9. AXIAL POST GAD LAVA BH (delay)**

OPTIONAL POST GAD LAVA AXIAL BH (10-15 min delay) (if suspect CholangioCa , PSC or slow-filling hemangioma)