

MRCP with Contrast

(revised 12/2009)

Maximum slice thickness 6mm

Maximum gap 1.5

- 1. SCOUT**
- 2. AXIAL T1 IN & OUT PHASE**
- 3. AXIAL TSE T2 resp comp (top of liver to bottom of pancreas)**
- 4. Thin 3mm CORONAL SSFSE T2 BH (through distal duct/pancreatic head)**
- 5. Thin 3mm AXIAL SSFSE T2 BH (above porta hepatis to bottom of pancreatic head)**
- 6. Thick SLAB MRCP – (12 slab 10°)**
- 7. 3D ASSET TSE navigated MRCP**
- 8 AXIAL T2 TSE FS resp comp (Cover top of liver to bottom of pancreas)**
- 9. Diffusion B500/B0**
- 10. AXIAL LAVA BH (Cover top of liver to bottom of pancreas)**

Inject contrast

- 11. AXIAL LAVA BH (multiple phases same as liver protocol) (Cover top of liver to bottom of pancreas)**
- 12. CORONAL LAVA (cover abdomen)**
- 13. AXIAL LAVA BH 5 minute delay**