

Osteomyelitis Lower Extremity

(revised 12/2008)

Contrast dosage 0.1mmol/kg

Mark area(s) of interest(s). Do not depress capsule into skin, place superficially

If area of interest ambiguous call/page the resident/attending

If poor T2 fat saturation use STIR

- 1. SCOUT**
- 2. SAG T1 SE**
- 3. COR T1 SE**
- 4. SAG STIR**
- 5. COR STIR**
- 6. AX T1**
- 7. AX T2 FS (TE 50-60)**
- 8. AX PRE T1 FAT SAT**

Inject Contrast

- 9. SAG POST T1 FAT SAT**
- 10. COR POST T1 FAT SAT**
- 11. AXIAL POST T1 FAT SAT**