

Penis

(revised March 22, 2012)

Use Large Loop Coil and Body array

Use small fields of view (~12-14cm)

Contrast dosage 0.1mmol/kg

Position penis parallel to sagittal plane, directed cephalically. Tape to abdomen. Elevate scrotum with towel.

Charges: MR Pevlis w/wo, contrast dose

1. SCOUT

2. AXIAL T2 SSFSE 6mm (Cover entire pelvis)

Small Field of View Penis

3. CORONAL TSE T2

4. SAGITTAL TSE T2

5 AXIAL T1 TSE FS

INJECT

6. AXIAL T1 TSE FS

7. CORONAL T1 TSE FS