

# Perirectal Fistula

(revised 3/2009)

Use body phased array coil.

Have patient empty bladder before scanning

Contrast dosage 0.1mmol/kg

**Max slice thickness 6mm, gap 1.5**

- 1. SCOUT BH**
- 2. AXIAL T2 (Large FOV entire pelvis)**
- 3. CORONAL SSFSE T2 NON BH  
(Large FOV entire pelvis)**
- 4. Thin AXIAL T2 TSE FS 3mm high  
resolution SMALL FOV (only rectal and  
anal area)**
- 5. PRE GAD AXIAL LAVA T1**

## **Inject contrast**

- 6. POST GAD AXIAL LAVA T1**
- 7. POST GAD COR LAVA T1**
- 8. POST GAD SAG LAVA T1**