

Renal

(revised 12/2009)

Body phase array coil.

Use power injector

Contrast dosage 0.1mmol/kg

Max slice thickness 6mm, gap 1.5

- 1. SCOUT BH**
- 2. CORONAL FSE T2 NON BH (cover kidneys)**
- 3. AXIAL FSE T2 NON BH (cover kidneys)**
- 4. AXIAL FSE T2 fs NON BH (cover kidneys)**
- 5. AXIAL IN and OPPOSED (cover kidneys)**
- 6. AX DIFFUSION B500, B0 (cover kidneys)**
- 7. CORONAL LAVA T1 BH (cover kidneys)**
- 8. AXIAL LAVA T1 BH (cover liver and kidneys)**

Inject contrast

- 9. AXIAL DYNAMIC POST GAD LAVA (cover liver and kidneys)**
- 10. SAGITTAL POST GAD LAVA (90 sec nephrograph kidney to kidney)**
- 11. CORONAL LAVA T1 BH (3 min excretory phase)**
- 13. AXIAL POST GAD CORONAL LAVA (4 min)**

Optional:

- CORONAL LAVA BH (6 min MR IVP phase)**
- SAGITTAL FSE or SSFSE T2**