

Scrotum/Testes

Body phased array coil and large loop coil.

Use a small FOV except for sequence (2)

Contrast dosage 0.1mmol/kg

1. SCOUT

2. AXIAL TSE T2 (Cover entire pelvis)

3. CORONAL TSE T2 (small fov)

4. SAGITTAL TSE T2 (small fov)

5. SAG TSE T1 (small fov)

6. AXIAL fs SPGR (small fov)

Inject Contrast

7. AXIAL fs T1 (small fov)

8. CORONAL fs T1 (small fov)

9. SAGITTAL fs T1 (small fov)