

# Urethral diverticulum

(revised March 22, 2012)

Use body phased array coil.

Have patient empty bladder before scanning

Contrast dosage 0.1mmol/kg

**Max total slice thickness 6mm, gap 1.5**

- 1. SCOUT BH**
- 2. CORONAL in and out of phase (Large FOV entire pelvis)**
- 3. AXIAL T2 fs (Large FOV entire pelvis)**
- 4. SAGITTAL T2 TSE SMALL FOV include uterus. pelvic sidewall to sidewall**
- 5. Thin AXIAL T2 TSE 3mm high resolution SMALL FOV (only urethral area – don't include uterus)**
- 6. PRE GAD AXIAL LAVA – Large FOV**

## Inject contrast

- 7. POST GAD AXIAL LAVA T1 – Large**
- 8. POST GAD COR LAVA T1 – Large**
- 9. POST GAD SAG LAVA T1 – Large**
- 10 POST GAD AXIAL T1 fs SMALL FOV thin slice 3mm urethral area only**