

Liver Donor

(revised 5/2018)

Body phase array coil.

Use power injector

Cover entire liver

Eovist Contrast

Siemens

Red = new

- 1.SCOUT BH**
- 2.COR HASTE**
- 3.AXIAL HASTE**
- 4.MRCP 12 thick slab**
- 5.MRCP 3D triggered**
- 6.AXIAL Dixon OUT/IN PHASE w liverlab**
- 7.AXIAL VIBE pre**

Inject Eovist– Fluoro-prep

- 8. AXIAL POST GAD VIBE BH dynamic**
- 9. COR POST GAD VIBE BH (2 min)**
- 10. AX POST GAD VIBE BH (delay)**
- 11. COR HASTE 3mm thin small FOV**
- 12. AX T2 fs trig blade**
- 13. AX TRUFISP**
- 14. DWI B50/500**
- 15. Ax VIBE 20 min**
- 16. Cor VIBE 20 min**

Call to check after MRCPs