

MRA/MRV Ablavar for

- AAA/dissection**
- mesenteric ischemia/SMA**
- renal artery**
- IVC**
- pelvic congestion syndrome**

(revised 3/2011, 12/2011)

Please do two station for AXIAL if large coverage

- 1. COR T2 SSFSE BH**
- 2. AX T1 IN/OUT PHASE BH**
- 3. AX T2 FIESTA BH**
- 4. Pre AX LAVA Low Res (5-6mm)**
- 5. Pre COR MRA MASK**

Inject Ablavar 0.03 mmol/kg, @1cc/sec, 30 cc saline

- 6. Post COR MRA**
- 7. Post AX LAVA Low Res**
- 8. Repeat Post COR MRA at 2min, 3min**
- 9. Post AX LAVA High Res (1-2mm) (after 5min)**
- 10. Post COR LAVA High Res (1-2mm)**

- 11. Optional TRICKS with 2cc @1cc/sec**